

## CROWN UNDER EXISTING PARTIAL



1. Prepare tooth using normal procedures. Retract tissue as indicated.



2. Seat partial fully to place over prepared tooth.

3. Using a flexible, light bodied impression or bite registration material, inject a "blob" all around clasp and prepared tooth, making

sure that it extends to the marginal area and fully encompasses all involved clasp arms and rests.



4. Remove partial; the "blob" will come out with it. Remove the "blob" and set aside.



5. Without partial in place and using standard procedure, take final impression for crown.

6. Send impression, "blob", opposing model, and bite registration to laboratory.



### Laboratory views

