


# BruxZir<sup>®</sup> is more brawn than beauty

You'll be impressed by the esthetics of BruxZir Solid Zirconia when prescribed instead of metal occlusal PFMs and full-cast metal restorations. Virtually chip-proof, BruxZir is a full-contour zirconia crown, bridge, inlay or onlay with no porcelain overlay.

Designed and milled using CAD/CAM technology, BruxZir restorations are sintered for 6.5 hours at 1,530 degrees Celsius. The final BruxZir crown or bridge emerges nearly "bulletproof" and is glazed to a smooth surface.


## Wear Compatibility



In a recent study at the University of Alabama at Birmingham to measure the volume metric loss of enamel, glazed BruxZir was found to be wear compatible with enamel and virtually identical to glazed IPS e.max.

IPS e.max is a registered trademark of Ivoclar Vivadent.

Polish BruxZir restorations with the porcelain polishing system of your choice. Or, purchase the Axis Dental BruxZir Adjustment & Polishing Set (LS-7579), developed by Dr. Michael DiTolla, through your local dental dealer or by calling 800-355-5063.



Before



BruxZir Solid Zirconia Crown



Full-Cast Gold Crown



Porcelain and Metal Crown